Medical Device Innovation through Biodesign Process

Adarsh Somayaji ED15B001

Company: Shira MedTech Pvt. Ltd.

Mentors: Anand Parikh, Prasann Jain



- To identify problems in a clinical environment that can be solved through technology
- To research on & validate the scope of the problems
- To shortlist problems & ideate on solution concepts for the same
- To develop a proof of concept prototype for the best solution

What is Biodesign?

Biodesign is an innovation model that begins with careful identification of a clinical need and moves in a stepwise approach through inventing & planning the implementation of a marketable solution

Identify

Invent

Implement

Strategic Focus

Clinical Immersion

Need Statement

Disease State Study

Needs Validation

Needs Filtering

Identify

Strategic Focus

The following departments were chosen as a part of our strategic focus:

- Emergency Medicine
- General Surgery
- Gastro-Intestinal Surgery
- Plastic Surgery
- Nephrology & Urosurgery
- Neurosurgery

Phase I - Dec

Phase II - Jan

Strategic Focus

Clinical Immersion

Need Statement

Disease State Study

Needs Validation

Needs Filtering

Identify

Clinical Immersion

Procedures Observed

~160

- All observations documented with time stamps
- Interactions with junior & senior doctors to understand the procedural pain points



Defective CBD stents extracted after an ERCP procedure

Strategic Focus

Clinical Immersion

Need Statement


Disease State Study

Needs Validation

Needs Filtering

Identify

Need Statement Development



Strategic Focus

Clinical Immersion

Need Statement


Disease State Study

Needs Validation

Needs Filtering

Identify

Disease State Study



Strategic Focus

Clinical Immersion

Need Statement

Disease State Study


Needs Validation

Needs Filtering

Identify



- Questionnaire circulated amongst doctors from different departments
- Doctors asked to rate each need statement based on 2 parameters
 - Patient/Provider Impact
 - Treatment Landscape



Strategic Focus

Clinical Immersion

Need Statement

Disease State Study


Needs Validation

Needs Filtering

Identify

Needs Filtering - Filter 1

The feedback scores from the questionnaire as we was used as a primary filter for the need statements




Market Analysis

Market Analysis involved analysis of the following for each need statement

- Epidemiology data to estimate market size & to segment the market
- Existing products to find gaps in user requirements & device specifications
- Opportunities & threats for developing solutions

Needs Filtering - Filter 2

Data from the market analysis & disease state studies was used to further filter the needs



Final Needs

A way to prevent blockage of **Endo-Tracheal** Tubes due to accumulation of secretions while reducing frequency of routine interventions

A way to provide prone ventilation to ARDS patients while accommodating for CPR A way to prevent needle prick related failures of Arteriovenous Fistulas for patients on dialysis

Final Needs

A way to prevent blockage of Endo-Tracheal Tubes due to accumulation of secretions while reducing frequency of routine interventions


A way to provide prone ventilation to ARDS patients while accommodating for CPR

A way to prevent needle prick related failures of Arteriovenous Fistulas



- ARDS: Acute Respiratory Distress Syndrome
- Inflammation of lung tissue leading to collapse of air sacs & accumulation of fluid in the lung
- Triggered by a variety of factors such as sepsis, bacterial or viral lung infections, inhalation of vomit or water, pancreatitis, pneumonia etc.


ARDS: The Numbers






- All ARDS patients are on mechanical ventilation while doctors treat root cause of ARDS
- Two methods of mechanical ventilation
 - Supine Ventilation
 - Prone Ventilation
 - Better for the patient as it is shown to result in improved oxygenation & lower chances of ventilator induced lung injury

Prone Ventilation: The Numbers



Prone Ventilation Physiology



Shape matching and gravity

2 effects balance each other in prone position

- Gravity
- Shape Matching

Shape matching

24



- Needs 6-8 trained personnel
- Difficult to provide CPR
- Detachment of tubings & lines
- Physical trauma to the patient while performing the procedure
- Difficult to examine patient in prone position
- Difficult to feed & suction patient
- Pressure sores & corneal Injuries

Need Criteria for MVP

Need Criteria - Must Haves
Prevention of physical trauma to patient
Provision for CPR
Accommodation of tracheostomy tubes
Prevention of dislodgement of attachments

Need Criteria for MVP

Need Criteria - Good to Haves
Reduction in effort in performing procedure
Reduction in time taken in performing procedure
Cost Effectiveness
Prevention of pressure sores
Prevention of Corneal injuries

I: Concept Generation

Invent

II: Concept Selection Brainstorming

Concept Screening

Concept Validation

Prototyping



Ideation & Brainstorming

Brainstorming Rules:

- Defer Judgement
- Encourage Wild Ideas
- Build on Ideas of Others
- Go for Quantity
- One Conversation at a Time
- Stay Focused
- Be Visual



3 sessions

132 ideas

I: Concept Generation Brainstorming

Concept Screening

Concept Validation

Prototyping

Invent

II: Concept Selection

Concept Screening

- Grouping & Organization of Ideas
- Filtering out infeasible concepts
- Concept Mapping

Concept Groups for flipping patient

Concept Groups for preventing trauma

Concept Groups for constraining tubings

I: Concept Generation Г

Brainstorming

Concept Screening

Concept Validation

Prototyping

Invent

II: Concept Selection



- Toy models made of lego, clay etc were developed for each main concept
- Models shown to ~15 intensivists at different hospitals for feedback
- New constraints of hospital economics found & applied to shortlist the best concept

I: Concept Generation Brainstorming

Concept Screening

Concept Validation

Prototyping

Invent

II: Concept Selection

References

- 1. Zenios, S., Makower, J., Yock, P., "Biodesign: The process of innovating medical technologies", Cambridge University Press 2010
- 2. Fan, E., Brodie, D., Slutsky, A.S., "Acute Respiratory Distress Syndrome: Advances in Diagnosis and Treatment", JAMA. 2018; 319(7):698–710
- 3. Mancebo, J., Fernández, R., Blanch, L., et al. "A Multicenter Trial of Prolonged Prone Ventilation in Severe Acute Respiratory Distress Syndrome", American Journal of Respiratory and Critical Care Medicine, 173(11), 1233–1239
- 4. Gattinoni, L., Taccone, P., Carlesso, E., and Marini, J.J., "Prone Position in Acute Respiratory Distress Syndrome. Rationale, Indications, and Limits" American Journal of Respiratory and Critical Care Medicine, 188(11)
- 5. http://ardsglobal.org/acute-respiratory-distress-syndrome/

Thank You!